

## Differential biochemical tests to characterize a species of bacteria

For every biochemical process they perform, they need to have the right enzyme(s)  
Enzyme expression is **species-specific**

Therefore, by looking at a **pattern** of biochemical activity, you can identify the organism.

Goal: determine biochemical profile of organism. No single test can tell you!

- Substrate use
  - End product formation
- 


### #13: Fermentation of sugars

Differential tests. Q: can UNK ferment the sugar? Choices: glucose, sucrose, lactose

Fermentation: ATP production without using O<sub>2</sub>.

**Obligate aerobes:** can't ferment

**Facultative anaerobes/ Indifferent/ Obligate anaerobes:** can ferment various sugars depending on species



First, look for turbidity (cell growth). If positive, ask:

- **acidic end products** produced?

pH indicator: Brom Cresol Purple (BCP glucose/sucrose/lactose)

7.0 purple <5.4 yellow (+ test) >8.0 blue (products of cellular respiration)

- **Gas end products (CO<sub>2</sub>, H<sub>2</sub>) may** also be produced.

Durham tube: collects gas as a bubble. In these tests, the gas is H<sub>2</sub>. CO<sub>2</sub>, when produced, is water soluble. H<sub>2</sub> is not.

Sucrose and lactose are **disaccharides**. In order to ferment these sugars, bacteria must possess enzymes for cleavage, and isomerization to glycolytic substrates.

In these tests: one sugar per tube. Ask: acid? (yellow=fermented); gas?

---

### #14: Triple Sugar Iron (TSI): see textbook table 6.6 & Dr. Metcalf's photos

Differential test. Very useful for distinguishing enteric bacteria (mostly gram – bacilli)

Key points:

- All THREE sugars are in this media.
- There is lots of oxygen present on the slant surface; less inside; and very little at the butt

Questions to ask of TSI results:

1. Do the bacteria ferment? If so, which sugars?

All 3 sugars are present in TSI, but 10x more sucrose & lactose than glucose

In presence of abundant S&L, even when oxygen is available, bacteria will preferentially ferment (if they can) and turn off cellular respiration pathway.

- In the butt: little or no oxygen, so only fermentative metabolism can occur. Yellow butt = acid pH = glucose fermented. No color change = no glucose fermentation

- In the slant: fermentation of high concentration sucrose and/or lactose if possible (yellow, acid); otherwise oxidative metabolism occurs (red, alkaline)
2. Cysteine metabolism: production of H<sub>2</sub>S end product, appears as a black ppt
  3. Gas production: cracks in medium

pH indicator is phenol red. 7.3 = red <6.8 = yellow (acid= + fermentation)

Where does fermentation occur? **Natural oxygen gradient in an agar slant:** lots on surface, little at bottom of tube.

Interpretation is best if you can compare colors of multiple tubes, including uninoculated control.

**Results expressed as slant/butt** A=acid, K=alkaline; gas or no gas; H<sub>2</sub>S or no

Interpretation:

Butt: yellow = glucose fermented; red or unchanged = can't ferment glucose

Slant: yellow (A) = other sugar(s) fermented

red (K) = cannot ferment others, glucose used aerobically, runs out of glucose and aerobically metabolizes proteins/aminoacids producing alkaline products (NH<sub>3</sub>)

### #15: Catalase & oxidase tests

Catalase and oxidase are enzymes containing a distinctive prosthetic group, or chemical attachment, called an **iron porphyrin** group. Classic example is hemoglobin. Others are various proteins of the electron transport chain, called cytochromes. The electron transport system is central to aerobic respiration, carrying electrons from NAD/FAD to oxygen, while making ATP. Another iron porphyrin protein is catalase, which breaks down toxic H<sub>2</sub>O<sub>2</sub> into H<sub>2</sub>O and O<sub>2</sub>.

If a cell can make catalase, it can make iron porphyrin groups and is expected to make cytochromes. A positive test for catalase, therefore, indicates cells which can respire aerobically. Test is especially useful for differentiating Gram + cocci.

There are MANY cytochromes, however. Different bacteria produce different ones. The oxidase test specifically detects cytochrome C, which is actually rare in bacteria. *Pseudomonas* genus is positive.

**Catalase -** = ferment only, cannot respire aerobically (THIS DOES NOT MEAN anaerobe); usually very small colonies as fermentation is less efficient than aerobic respiration

**Catalase +** = can make iron porphyrin groups; can make cytochromes & respire aerobically

**Oxidase +** = all of above plus specifically makes cytochrome C

**Oxidase -** = no cytochrome c; ???catalase or aerobic respiration

### Inoculate:

#13: Using unk SLANT & loop (not broth as it says in lab manual), inoculate glucose, lactose & sucrose broths (make sure no gas bubble present). CHECK for bubbles in Durham tube (shouldn't be any); bring labels with you so you don't mix up the sugars; use aseptic technique

#14: Using unk slant with needle, inoculate TSI stab then streak. Aim for center of media.

#15: **Streak for colony isolation:** NAG plate with unk from BROTH